


Anti-Hsp20/Hspb6 Rabbit Polyclonal Antibody For Wb, Elisa - O14558

Item Number: CM0029088



Introduction

CamelBio supplies a rabbit polyclonal antibody targeting human HSP20/HSPB6 (Hsp20, HEL55), validated for WB and ELISA. Cross-reacts with mouse. Ideal for studying molecular chaperone function, cardiac contractility, and stress response pathways.

[Learn More](#)

Parameter	Value
Product Name	HSP20/HSPB6 Rabbit pAb
Remarks/Alias	HEL55; Hsp20; PPP1R91; HSP20/HSPB6
Species	Human
GeneID (Human)	126393
GeneID	126393
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human HSP20/HSPB6 (NP_653218.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O14558
Protein Weight	17kDa
Shipping	Ice bag