


## Anti-Hsp27/Hspb1 Monoclonal Antibody For Wb - P04792

Item Number: CM0003964



### Introduction

Rabbit monoclonal antibody targeting HSP27/HSPB1 (P04792). Validated for WB, IF/ICC, IF-P, ELISA. Reacts with human, mouse, rat. Ideal for stress response, chaperone, and neurodegeneration research.

[Learn More](#)

| Parameter        | Value                                                                      |
|------------------|----------------------------------------------------------------------------|
| Product Name     | HSP27/HSPB1 Rabbit mAb                                                     |
| Remarks/Alias    | CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; HSP27/HSPB1 |
| Species          | Human                                                                      |
| GeneID (Human)   | 3315                                                                       |
| GeneID           | 3315                                                                       |
| Immunogen        | Synthetic Peptide                                                          |
| Source           | Rabbit                                                                     |
| Category         | Monoclonal Antibodies                                                      |
| Application      | WB, IF/ICC, IF-P, ELISA                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                          |
| SWISS            | P04792                                                                     |
| Protein Weight   | 23kDa                                                                      |
| Shipping         | Ice bag                                                                    |