


Anti-Hsp27/Hspb1 Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ip, If/Icc, Ihc-P, Elisa - P04792

Item Number: CM0009649



Introduction

KD validated rabbit monoclonal antibody against human HSP27/HSPB1, suitable for WB, IP, IF/ICC, IHC-P, and ELISA. Reliable tool for studying molecular chaperone function, stress response, and cytoskeleton dynamics.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] HSP27/HSPB1 Rabbit mAb
Remarks/Alias	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102
Species	Human
GeneID (human)	3315
GeneID	3315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human HSP27/HSPB1 (NP_001531.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P04792
Protein Weight	23kDa
Shipping	Ice bag