

Anti-Hsp40 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P25685

Item Number: CM0004240



Introduction

High-performance rabbit monoclonal antibody targeting human HSP40/DNAJB1 for WB, intracellular FC, and ELISA. Cross-reacts with human, mouse, rat. Ideal for chaperone, heat shock, and HSC70 interaction studies.

[Learn More](#)

Parameter	Value
Product Name	HSP40 Rabbit mAb
Remarks/Alias	DNAJB1; HSPF1; Hdj1; Hsp40; RSPH16B; Sis1; dnaJ homolog subfamily B member 1; HSP40
Species	Human
Gene ID (Human)	3337
Gene ID	3337
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 241-340 of human DNAJB1/HSP40 (NP_006136.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25685
Protein Weight	27kDa/38kDa
Shipping	Ice bag