

Anti-Hsp47/Serpinh1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50454

Item Number: CM0003359

Introduction

High-quality anti-HSP47/SERPINH1 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human and mouse. Targets serpin H1 collagen chaperone for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	HSP47/SERPINH1 Rabbit mAb
Remarks/Alias	CBP1; CBP2; OI10; gp46; AsTP3; HSP47; PIG14; PPROM; RA-A47; SERPINH2; HSP47/SERPINH1
Species	Human
GeneID (Human)	871
GeneID	871
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 319-418 of human HSP47/SERPINH1 (P50454).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50454
Protein Weight	46kDa
Shipping	Ice bag