

Anti-Hsp47/Serpinh1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P50454

Item Number: CM0014979

Introduction

Rabbit polyclonal antibody targeting human HSP47 (SERPINH1, P50454) for Western blot, IF/ICC, and ELISA. Detects endogenous HSP47 in human and mouse samples; ideal for collagen chaperone studies.

[Learn More](#)

Parameter	Value
Product Name	HSP47/SERPINH1 Rabbit pAb
Remarks/Alias	CBP1; CBP2; OI10; gp46; AsTP3; HSP47; PIG14; PPROM; RA-A47; SERPINH2; HSP47/SERPINH1
Species	Human
GeneID (Human)	871
GeneID	871
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 139-418 of human HSP47/SERPINH1 (NP_001226.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50454
Protein Weight	46kDa
Shipping	Ice bag