


Anti-Hsp60/Hspd1 Monoclonal Antibody (Rabbit Mab) For Wb, Ihc-P, If/Icc, Elisa - P10809

Item Number: CM0003113



Introduction

High-quality anti-HSP60/HSPD1 rabbit monoclonal antibody for WB, IHC-P, IF/ICC, and ELISA. Validated in human, mouse, rat, and wheat; targets mitochondrial chaperonin HSP60 (61 kDa) for protein folding research.

[Learn More](#)

Parameter	Value
Product Name	HSP60/HSPD1 Rabbit mAb
Remarks/Alias	HLD4; CPN60; GROEL; HSP60; HSP65; SPG13; HSP-60; HuCHA60; HSP60/HSPD1
Species	Human
GeneID (human)	3329
GeneID	3329
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Wheat
SWISS	P10809
Protein Weight	61kDa
Shipping	Ice bag