


Anti-Hsp60/Hspd1 Rabbit Polyclonal Antibody For Wb - P10809

Item Number: CM0015405



Introduction

High-quality rabbit polyclonal antibody targeting human HSP60/HSPD1. Validated for WB, IF/ICC, IHC-P, IP, ELISA; cross-reacts with mouse, rat. Suitable for mitochondrial chaperone research.

[Learn More](#)

Parameter	Value
Product Name	HSP60/HSPD1 Rabbit pAb
Remarks/Alias	HLD4; CPN60; GROEL; HSP60; HSP65; SPG13; HSP-60; HuCHA60; HSP60/HSPD1
Species	Human
GeneID (Human)	3329
GeneID	3329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-240 of human HSP60/HSPD1 (NP_002147.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10809
Protein Weight	61kDa
Shipping	Ice bag