

Anti-Hsp70 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P0Dmv8 / P0Dmv9

Item Number: CM0006886



Introduction

Recombinant rabbit monoclonal antibody against human HSP70 (HSPA1A/HSPA1B). Suitable for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat, wheat, maize. Immunogen: amino acids 500-600 of HSP70.

[Learn More](#)

Specification	Value
Product Name	HSP70 Rabbit mAb
Remarks/Alias	HSP70; HSP72; HSPA1; HSP70I; HSP70-1; HSP70-2; HSP70.1; HSP70.2; HSP70-1A; HEL-S-103
Species	Human
GeneID	3303/3304
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-600 of human HSP70. (NP_005336.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Wheat, Maize
SWISS	P0DMV8/P0DMV9
Protein Weight	70 kDa
Shipping	Ice bag