

Anti-Hsp70 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P0Dmv8 / P0Dmv9

Item Number: CM0002761



Introduction

Recombinant rabbit monoclonal antibody targeting human HSP70 (HSPA1A/HSPA1B). Validated for Western blot, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for studying chaperone functions, stress response, and protein quality control pathways.

[Learn More](#)

Parameter	Value
Product Name	HSP70 Rabbit mAb
Remarks/Alias	HEL-S-103; HSP70-1; HSP70-1A; HSP70.1; HSP70I; HSP72; HSPA1; HSP70
Species	Human
GeneID (Human)	3303/3304
GeneID	3303/3304
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 500-600 of human HSP70 (NP_005336.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DMV8/P0DMV9
Protein Weight	70 kDa
Shipping	Ice bag