

Anti-Hsp70/Hspa1A Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P0Dmv8/P0Dmv9

Item Number: CM0014129



Introduction

Rabbit polyclonal antibody against human HSP70 (HSPA1A), validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Targets the key molecular chaperone involved in stress response and protein folding.

[Learn More](#)

Parameter	Value
Product Name	HSP70 Rabbit pAb
Remarks/Alias	HSPA1A; HEL-S-103; HSP70-1; HSP70-1A; HSP70.1; HSP70I; HSP72; HSPA1; HSP70
Species	Human
GeneID (Human)	3303/3304
Immunogen	Synthetic peptide: a synthetic peptide corresponding to a sequence within amino acids 500-600 of human HSP70/HSPA1 (NP_005336.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DMV8/P0DMV9
Protein Weight	70 kDa
Shipping	Ice bag