


Anti-Hsp90A Polyclonal Antibody For Wb, If/Icc, Elisa - P07900

Item Number: CM0013987



Introduction

Rabbit polyclonal Hsp90 α antibody validated for Western blot, IF/ICC, and ELISA in human, mouse, and rat samples. Targets heat shock protein HSP 90-alpha, a molecular chaperone involved in cell cycle control and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	Hsp90 α Rabbit pAb
Remarks/Alias	EL52; HSPN; LAP2; HSP86; HSPC1; HSPCA; Hsp89; Hsp90; LAP-2; HSP89A; HSP90A; HSP90N; Hsp103; HSPCAL1; HSPCAL4; HEL-S-65p; Hsp90 α
Species	Human
GeneID (human)	3320
GeneID	3320
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07900
Protein Weight	85kDa
Shipping	Ice bag