


Anti-Hsp90A/B Rabbit Monoclonal Antibody For Wb, Elisa - P07900 / P08238

Item Number: CM0010290



Introduction

High-quality rabbit monoclonal antibody against human Hsp90 α / β , validated for WB and ELISA. Suitable for chaperone research, signal transduction, and mitochondrial import studies. Cross-reacts with human, mouse, rat, monkey.

[Learn More](#)

Parameter	Value
Product Name	Hsp90 α / β Rabbit mAb
Remarks/Alias	EL52; HEL-S-65p; HSP86; HSP89A; HSP90A; HSP90N; HSPC1; HSPCA; HSPCAL1; HSPCAL4; HSPN; Hsp103; Hsp89; Hsp90; LAP-2; LAP2; Hsp90 α / β
Species	Human
Gene ID (Human)	3320/3326
Gene ID	3320/3326
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 533-732 of human Hsp90 α (P07900)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P07900/P08238
Protein Weight	90kDa
Shipping	Ice bag