


Anti-Hsp90B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08238

Item Number: CM0024163



Introduction

Hsp90β Rabbit pAb validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human HSP90-beta (P08238), cross-reacts with mouse and rat. Ideal for studying chaperone-mediated signaling and stress responses.

[Learn More](#)

Parameter	Specification
Product Name	Hsp90β Rabbit pAb
Remarks/Alias	HSP84; HSPC2; HSPCB; D6S182; HSP90B; Hsp90β
Species	Human
GeneID (Human)	3326
GeneID	3326
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Hsp90β (NP_031381.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08238
Protein Weight	83kDa
Shipping	Ice bag