


Anti-Hsp90B Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08238

Item Number: CM0006425



Introduction

Hsp90β Rabbit monoclonal antibody validated for Western blot, immunohistochemistry, immunofluorescence, and ELISA. Cross-reacts with human, mouse, rat, and rice. Ideal for heat shock protein and chaperone studies.

[Learn More](#)

Parameter	Value
Product Name	Hsp90β Rabbit mAb
Remarks/Alias	HSP84; HSPC2; HSPCB; D6S182; HSP90B; Hsp90β
Species	Human
GeneID (Human)	3326
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Hsp90β. (NP_031381.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Rice
SWISS	P08238
Protein Weight	83kDa
Shipping	Ice bag