


Anti-Hsp90B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P08238

Item Number: CM0018086



Introduction

Anti-Hsp90 β Rabbit Polyclonal Antibody (KO Validated) for Western blotting, IHC-P, and ELISA; reactive with human, mouse, rat. Essential tool for studying molecular chaperone function in signal transduction, stress response, and cell cycle control.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Hsp90 β Rabbit pAb
Remarks/Alias	HSP84; HSPC2; HSPCB; D6S182; HSP90B; 0 β
Species	Human
GeneID (Human)	3326
GeneID	3326
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Hsp90 β (NP_031381.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08238
Protein Weight	83kDa
Shipping	Ice bag