

Anti-Hspa12A Monoclonal Antibody For Wb And Elisa - O43301

Item Number: CM0003652

Introduction

Rabbit monoclonal antibody against human HSPA12A, suitable for Western blot and ELISA. Detects endogenous HSPA12A in human, mouse, and rat samples. Ideal for studying SORL1-mediated trafficking and heat shock response.

[Learn More](#)

Parameter	Value
Product Name	HSPA12A Rabbit mAb
Remarks/Alias	HSPA12A ; HSPA12A
Species	Human
GeneID (Human)	259217
GeneID	259217
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human HSPA12A (O43301).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43301
Protein Weight	75 kDa
Shipping	Ice bag