

Anti-Hspa14 Rabbit Polyclonal Antibody For Wb, Elisa - Q0Vdf9

Item Number: CM0029354

Introduction

Rabbit polyclonal antibody against human HSPA14, validated for WB and ELISA; cross-reacts with mouse. For ribosome-associated complex and protein folding studies. Order from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	HSPA14 Rabbit pAb
Remarks/Alias	HSP70-4; HSP70L1; MSANTD7; HSPA14
Species	Human
Gene ID (Human)	51182
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-509 of human HSPA14 (NP_057383.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q0VDF9
Protein Weight	55kDa
Shipping	Ice bag