

Anti-Hspa1L Polyclonal Antibody For Wb - P34931

Item Number: CM0018908



Introduction

Anti-HSPA1L rabbit polyclonal antibody for WB and ELISA, recognizing human and mouse HSPA1L. Ideal for studying molecular chaperone functions, protein quality control, and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	HSPA1L Rabbit pAb
Remarks/Alias	HSP70T; hum70t; HSP70-1L; HSP70-HOM; HSPA1L
Species	Human
GeneID (Human)	3305
GeneID	3305
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 342-641 of human HSPA1L (NP_005518.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P34931
Protein Weight	70kDa
Shipping	Ice bag