


Anti-Hspa2 Monoclonal Antibody For Wb, If/Icc, Elisa - P54652

Item Number: CM0012683



Introduction

High-quality rabbit monoclonal antibody targeting human HSPA2 (HSP70-2), suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for investigating chaperone-mediated protein quality control and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	HSPA2 Rabbit mAb
Remarks/Alias	HSP70-2; HSP70-3; HSPA2
Species	Human
GeneID (Human)	3306
GeneID	3306
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54652
Protein Weight	70kDa
Shipping	Ice bag