


Anti-Hspa2 Rabbit Polyclonal Antibody - P54652

Item Number: CM0028942



Introduction

HSPA2 Rabbit pAb is a polyclonal antibody sourced from rabbit, validated for Western Blot (WB), IF/ICC, and ELISA. It reacts with human, mouse, and rat HSPA2, a heat shock protein involved in chaperone functions, protein quality control, and spermatogenesis. Ideal for stress response and cell biology studies.

[Learn More](#)

Parameter	Value
Product Name	HSPA2 Rabbit pAb
Remarks/Alias	HSP70-2, HSP70-3, HSPA2
Species	Human
GeneID (Human)	3306
GeneID	3306
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-639 of human HSPA2 (NP_068814.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54652
Protein Weight	70kDa
Shipping	Ice bag