

Anti-Hspa4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P34932

Item Number: CM0002807

Introduction

High-specificity rabbit monoclonal antibody against human HSPA4 (heat shock 70 kDa protein 4). Validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, rat. Suitable for stress response and protein folding studies.

[Learn More](#)

Parameter	Value
Product Name	HSPA4 Rabbit mAb
Remarks/Alias	RY; APG-2; HSPH2; hsp70; hsp70RY; HEL-S-5a; HS24/P52; HSPA4
Species	Human
GeneID (Human)	3308
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 700-800 of human HSPA4 (P34932).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34932
Protein Weight	94 kDa
Shipping	Ice bag