

Anti-Hspa4 Rabbit Polyclonal Antibody For Wb, Elisa - P34932

Item Number: CM0023863

Introduction

Rabbit polyclonal antibody against human HSPA4 (heat shock 70 kDa protein 4, APG-2). Validated for Western blot and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for stress biology and chaperone research.

[Learn More](#)

Parameter	Value
Product Name	HSPA4 Rabbit pAb
Remarks/Alias	RY; APG-2; HSPH2; hsp70; hsp70RY; HEL-S-5a; HS24/P52; HSPA4
Species	Human
GeneID (Human)	3308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 680-840 of human HSPA4 (NP_002145.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34932
Protein Weight	94kDa
Shipping	Ice bag