

Anti-Hspa4L Polyclonal Antibody For Wb, If/Icc, Elisa - 095757

Item Number: CM0021622

Introduction

High-quality HSPA4L rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, rat. Ideal for studies on heat shock protein 70 family chaperone function and stress response. Research use only.

[Learn More](#)

Parameter	Value
Product Name	HSPA4L Rabbit pAb
Remarks/Alias	APG1; APG-1; HSPH3; Osp94; HSPA4L
Species	Human
GeneID (Human)	22824
GeneID	22824
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 500-600 of human HSPA4L (NP_055093.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95757
Protein Weight	95kDa
Shipping	Ice bag