

Anti-Hspb2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16082

Item Number: CM0024063



Introduction

Anti-HSPB2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Immunogen: recombinant protein (aa 1-182). Suitable for stress response and DMPK regulation research.

[Learn More](#)

Parameter	Value
Product Name	HSPB2 Rabbit pAb
Remarks/Alias	MKBP, HSP27, Hs:78846, LOH11CR1K, HSPB2
Species	Human
GeneID (Human)	3316
GeneID	3316
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of human HSPB2 (NP_001532.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16082
Protein Weight	20kDa
Shipping	Ice bag