

Anti-Hspb7 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uby9

Item Number: CM0021564



Introduction

Rabbit polyclonal antibody targeting human HSPB7 (cvHSP), suitable for Western blot, IHC-P, and ELISA applications. Cross-reacts with mouse and rat. Supports cardiovascular and stress response research.

[Learn More](#)

Parameter	Value
Product Name	HSPB7 Rabbit pAb
Remarks/Alias	cvHSP; HSPB7
Species	Human
Gene ID (Human)	27129
Gene ID	27129
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human HSPB7 (NP_055239.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBY9
Protein Weight	19kDa
Shipping	Ice bag