

Anti-Hspb8/Hsp22 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Ujy1

Item Number: CM0002046



Introduction

High-quality rabbit monoclonal antibody against human HSPB8/HSP22, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for chaperone-assisted autophagy and muscle disease research.

[Learn More](#)

Parameter	Value
Product Name	HSPB8/HSP22 Rabbit mAb
Remarks/Alias	H11; HMN2; CMT2L; DHMN2; E2IG1; HMN2A; HSP22; CRYAC; HSPB8/HSP22
Species	Human
GeneID (Human)	26353
GeneID	26353
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HSPB8/HSP22 (Q9UJY1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJY1
Protein Weight	22kDa
Shipping	Ice bag