

Anti-Hspb8/Hsp22 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9Ujy1

Item Number: CM0011577



Introduction

Recombinant anti-HSPB8/HSP22 rabbit monoclonal antibody validated for WB, IF/ICC, IHC-P, ELISA. Targets human HSPB8; cross-reacts with mouse, rat, monkey. Suitable for protein quality control research.

[Learn More](#)

Parameter	Value
Product Name	HSPB8/HSP22 Rabbit mAb
Remarks/Alias	H11; HMN2; CMT2L; DHMN2; E2IG1; HMN2A; HSP22
Species	Human
GeneID (Human)	26353
GeneID	26353
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-196 of human HSPB8/HSP22 (NP_055180.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	Q9UJY1
Protein Weight	22kDa
Shipping	Ice bag