

Anti-Hspb8/Hsp22 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Ujy1

Item Number: CM0020591



Introduction

Rabbit polyclonal antibody against human HSPB8/HSP22 (H11, HMN2, CMT2L) validated in WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying protein quality control, muscle biology, and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	HSPB8/HSP22 Rabbit pAb
Remarks/Alias	H11; HMN2; CMT2L; DHMN2; E2IG1; HMN2A; HSP22; HSPB8/HSP22
Species	Human
GeneID (Human)	26353
GeneID	26353
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-196 of human HSPB8/HSP22 (NP_055180.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJY1
Protein Weight	22kDa
Shipping	Ice bag