


Anti-Hspe1/Hsp10/Cpn10 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P61604

Item Number: CM0011601



Introduction

Recombinant rabbit monoclonal antibody against human HSPE1 (HSP10/CPN10) for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for mitochondrial chaperone and stress response research.

[Learn More](#)

Parameter	Value
Product Name	HSPE1/HSP10/CPN10 Rabbit mAb
Remarks/Alias	EPF; CPN10; GROES; HSP10; HSPE1/HSP10/CPN10
Species	Human
Gene ID (Human)	3336
Gene ID	3336
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-150 of human HSPE1/HSP10/CPN10 (P61604)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61604
Protein Weight	11kDa
Shipping	Ice bag