

Anti-Htr1B Rabbit Polyclonal Antibody For Wb, Elisa - P28222

Item Number: CM0019137



Introduction

Rabbit polyclonal antibody against human HTR1B (5-HT1B serotonin receptor), validated for Western blot and ELISA applications with cross-reactivity in mouse and rat. Ideal for neuroscience and receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	HTR1B Rabbit pAb
Remarks/Alias	S12, 5-HT1B, HTR1D2, HTR1DB, 5-HT-1B, 5-HT1DB, 5-HT-1D-beta, HTR1B
Species	Human
Gene ID (Human)	3351
Gene ID	3351
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-390 of human HTR1B (P28222).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28222
Protein Weight	44kDa
Shipping	Ice bag