


Anti-Htr1D Rabbit Polyclonal Antibody For Wb, Elisa - P28221

Item Number: CM0020528



Introduction

High-quality rabbit polyclonal antibody targeting human HTR1D (serotonin receptor 1D/5-HT1D). Validated for WB and ELISA applications, with cross-reactivity to mouse and rat. Ideal for neuroscience and GPCR signaling studies.

[Learn More](#)

Parameter	Value
Product Name	HTR1D Rabbit pAb
Remarks/Alias	HTRL; RDC4; HT1DA; 5-HT1D; HTR1DA; HTR1D
Species	Human
GeneID (Human)	3352
GeneID	3352
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human HTR1D (NP_000855.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28221
Protein Weight	42kDa
Shipping	Ice bag