


Anti-Htr2A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P28223

Item Number: CM0023692



Introduction

Rabbit polyclonal antibody targeting human HTR2A (5-HT2A serotonin receptor). Validated for WB, IF/ICC, and ELISA in human and mouse samples. Ideal for neuroscience and signaling research.

[Learn More](#)

Parameter	Value
Product Name	HTR2A Rabbit pAb
Remarks/Alias	HTR2; 5-HT2A; HTR2A
Species	Human
GeneID (Human)	3356
GeneID	3356
Immunogen	Synthetic peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HTR2A (NP_000612.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P28223
Protein Weight	53kDa
Shipping	Ice bag