


Anti-Htr2B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P41595

Item Number: CM0026316



Introduction

High-quality rabbit polyclonal antibody against human serotonin receptor 2B (HTR2B), validated for WB, IF/ICC, and ELISA. Ideal for studying GPCR signaling, neurobiology, and cardiovascular development.

[Learn More](#)

Parameter	Specification
Product Name	HTR2B Rabbit pAb
Remarks/Alias	5-HT2B; 5-HT-2B; 5-HT(2B); HTR2B
Species	Human
Gene ID (Human)	3357
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 382-481 of human HTR2B (NP_000858.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41595
Protein Weight	54 kDa
Shipping	Ice bag