

Anti-Htr2C (5-Ht2C) Polyclonal Antibody For Wb, Elisa - P28335

Item Number: CM0016486



Introduction

Rabbit polyclonal antibody against human HTR2C (5-HT2C serotonin receptor), validated for Western blot and ELISA. Reacts with human, mouse, and rat; immunogen located in amino acids 120-220. Ideal for neuropharmacology research.

[Learn More](#)

Parameter	Value
Product Name	HTR2C Rabbit pAb
Remarks/Alias	HTR1C; 5-HT1C; 5-HT2C; 5HTR2C; 5-HTR2C; HTR2C
Species	Human
GeneID (human)	3358
GeneID	3358
Immunogen	Synthetic peptide within amino acids 120-220 of human HTR2C (NP_000859.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28335
Protein Weight	52kDa
Shipping	Ice bag