


Anti-Htr3A Monoclonal Antibody For Wb, If, Ihc-P, Elisa - P46098

Item Number: CM0007026



Introduction

Rabbit monoclonal anti-HTR3A antibody validated for western blot, immunofluorescence, immunohistochemistry, and ELISA. Detects serotonin receptor 3A in human, mouse, and rat samples; ideal for neuroscience and ion channel studies.

[Learn More](#)

Parameter	Value
Product Name	HTR3A Rabbit mAb
Remarks/Alias	HTR3; 5HT3R; 5-HT-3; 5-HT3A; 5-HT3R; HTR3A
Species	Human
GeneID (Human)	3359
GeneID	3359
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46098
Protein Weight	55kDa
Shipping	Ice bag