


Anti-Htr3A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P46098

Item Number: CM0030139



Introduction

Rabbit polyclonal antibody targeting human HTR3A (5-HT3 receptor), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying serotonin-gated ion channels in neuronal signaling.

[Learn More](#)

Parameter	Value
Product Name	HTR3A Rabbit pAb
Remarks/Alias	HTR3; 5HT3R; 5-HT-3; 5-HT3A; 5-HT3R; HTR3A
Species	Human
Gene ID (Human)	3359
Gene ID	3359
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46098
Protein Weight	55kDa
Shipping	Ice bag