

Anti-Htr4 Rabbit Polyclonal Antibody For Wb, Elisa - Q13639

Item Number: CM0025207



Introduction

High-quality rabbit polyclonal antibody targeting the human serotonin receptor 4 (HTR4). Validated for Western blotting and ELISA, with cross-reactivity to mouse and rat. Ideal for neuroscience, cardiovascular, and GPCR signaling studies.

[Learn More](#)

Parameter	Value
Product Name	HTR4 Rabbit pAb
Remarks/Alias	5-HT4; 5-HT4R; HTR4
Species	Human
GeneID (Human)	3360
GeneID	3360
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-384 of human HTR4 (NP_000861.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13639
Protein Weight	44kDa
Shipping	Ice bag