


Anti-Htra2 Rabbit Monoclonal Antibody For Wb, Elisa - O43464

Item Number: CM0011282



Introduction

High-quality anti-HTRA2 rabbit monoclonal antibody validated for Western blot (WB) and ELISA. Recognizes human, mouse, and rat HTRA2 (~49 kDa). Ideal for apoptosis, neurodegeneration, and Parkinson disease research.

[Learn More](#)

Parameter	Value
Product Name	HTRA2 Rabbit mAb
Remarks/Alias	OMI; MGCA8; PARK13; PRSS25; HTRA2
Species	Human
Gene ID	27429
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HTRA2 (O43464)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43464
Protein Weight	49kDa
Shipping	Ice bag