


Anti-Htra2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - O43464

Item Number: CM0014230



Introduction

Rabbit polyclonal antibody targeting human HTRA2 (OMI/PARK13) serine protease, suitable for WB, IHC-P, IP, and ELISA applications. Cross-reacts with mouse and rat. Useful for apoptosis and Parkinson's disease studies. Uniprot O43464.

[Learn More](#)

Parameter	Value
Product Name	HTRA2 Rabbit pAb
Remarks/Alias	OMI; MGCA8; PARK13; PRSS25; HTRA2
Species	Human
GeneID (Human)	27429
GeneID	27429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 334-458 of human HTRA2 (NP_037379.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43464
Protein Weight	49kDa
Shipping	Ice bag