

Anti-Human Agtr2 Monoclonal Antibody For Flow Cytometry - P50052

Item Number: CM0010704



Introduction

-conjugated rabbit monoclonal antibody targeting human AGTR2 (AT2 receptor) for flow cytometry. Suitable for cardiovascular and cell signaling research. Cross-reacts with human AGTR2.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human AGTR2 mAb
Remarks/Alias	AT2; ATGR2; MRX88; AGTR2
Species	Human
GeneID (Human)	186
GeneID	186
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-45 of human AGTR2 (NP_000677.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P50052
Protein Weight	41kDa
Shipping	Ice bag