


Anti-Human Albumin Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02768

Item Number: CM0018334



Introduction

CamelBio's anti-human Albumin rabbit polyclonal antibody is validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Supports research on albumin function, zinc transport, and drug binding.

[Learn More](#)

Parameter	Value
Product Name	Albumin Rabbit pAb
Remarks/Alias	HSA; FDAHT; PRO0883; PRO0903; PRO1341; Albumin
Species	Human
GeneID (human)	213
GeneID	213
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-609 of human Albumin (NP_000468.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02768
Protein Weight	69kDa
Shipping	Ice bag