

Anti-Human B7-H2/Cd275 Monoclonal Antibody For Fc - 075144

Item Number: CM0008633



Introduction

Rabbit monoclonal antibody against human B7-H2/CD275 (ICOSLG), validated for flow cytometry. Suitable for T-cell costimulation and immune checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human B7-H2/CD275 mAb
Remarks/Alias	ICOSLG; B7-H2; B7H2; B7RP-1; B7RP1; CD275; GL50; ICOS-L; ICOSL; LICOS; ICOS ligand
Species	Human
GeneID (Human)	23308
GeneID	23308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-256 of human B7-H2/CD275 (NP_056074.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75144
Protein Weight	20kDa/33kDa/34kDa
Shipping	Ice bag