


Anti-Human Cd100/Sema4D Monoclonal Antibody () For If/Icc And Fc - Q92854

Item Number: CM0010698



Introduction

-conjugated rabbit monoclonal antibody targeting human Semaphorin-4D (CD100/SEMA4D) for immunofluorescence (IF/ICC) and flow cytometry (FC) applications. This antibody recognizes human SEMA4D, a cell surface receptor regulating neuronal development and immune responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD100/SEMA4D mAb
Remarks/Alias	SEMA4D; C9orf164; CD100; M-sema-G; SEMAJ; coll-4; semaphorin-4D
Species	Human
GeneID (Human)	10507
GeneID	10507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-734 of human Semaphorin-4D/CD100 (NP_006369.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	Q92854
Protein Weight	82kDa/96kDa
Shipping	Ice bag