

Anti-Human Cd116 Monoclonal Antibody For Flow Cytometry - P15509

Item Number: CM0007733



Introduction

-conjugated rabbit monoclonal antibody targeting human CD116 (CSF2RA/GM-CSFR α). Validated for flow cytometry, suitable for studying hematopoietic cell proliferation and differentiation. SWISS: P15509.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD116 mAb
Remarks/Alias	GMR; CD116; CSF2R; SMDP4; CDw116; CSF2RX; CSF2RY; GMCSFR; CSF2RAX; CSF2RAY; alphaGMR; GMR-alpha; GMCSFR-alpha; GM-CSF-R-alpha
Species	Human
GeneID (Human)	1438
GeneID	1438
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-320 of human CD116 (NP_006131.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15509
Protein Weight	26kDa/31kDa/38kDa/43kDa/46kDa/49kDa
Shipping	Ice bag