

Anti-Human Cd117/C-Kit Monoclonal Antibody For Flow Cytometry - P10721

Item Number: CM0012391



Introduction

-conjugated rabbit monoclonal antibody targeting human CD117/c-Kit, validated for flow cytometry. Suitable for hematopoietic stem cell and mast cell research, and oncology applications.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD117/c-Kit mAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC
Species	Human
GeneID (Human)	3815
GeneID	3815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-520 of human CD117/c-Kit (NP_000213.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag