

Anti-Human Cd11B Monoclonal Antibody For Fc - P11215

Item Number: CM0008115



Introduction

Mouse monoclonal antibody against human CD11b (Integrin alpha-M, Mac-1) validated for flow cytometry. Ideal for studying complement receptor functions, neutrophil adhesion, and monocyte/macrophage biology. SWISS: P11215.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD11b mAb
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6
Species	Human
GeneID (human)	3684
GeneID	3684
Immunogen	Synthetic peptide (A synthesized peptide derived from human CD11b)
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P11215
Protein Weight	127kDa
Shipping	Ice bag