


Anti-Human Cd123 (IL3R α) Monoclonal Antibody For Flow Cytometry - P26951

Item Number: CM0009158



Introduction

-conjugated rabbit monoclonal antibody targeting human CD123, the IL-3 receptor alpha subunit. Validated for flow cytometry, suitable for studying hematopoietic progenitor cells and IL-3 signaling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD123 mAb
Remarks/Alias	IL3R; CD123; IL3RX; IL3RY; IL3RAY; hIL-3Ra; IL-3R-alpha
Species	Human
GeneID (Human)	3563
GeneID	3563
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26951
Protein Weight	43 kDa
Shipping	Ice bag