

Anti-Human Cd124/Il-4Rα Monoclonal Antibody For Fc - P24394

Item Number: CM0010477



Introduction

-conjugated recombinant rabbit monoclonal antibody targeting human CD124/IL-4Rα. For flow cytometry research on IL-4/IL-13 signaling, Th2 differentiation, and allergic inflammation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD124/IL-4Rα mAb
Remarks/Alias	CD124; IL4RA; IL-4RA
Species	Human
Gene ID (Human)	3566
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-232 of human CD124/IL-4Rα (NP_000409.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P24394
Protein Weight	26kDa/88kDa/90kDa
Shipping	Ice bag