

Anti-Human Cd124/IL-4R α Monoclonal Antibody For Flow Cytometry - P24394

Item Number: CM0011802



Introduction

-conjugated rabbit monoclonal antibody targeting human CD124/IL-4R α , validated for flow cytometry. Suitable for detecting the IL-4 receptor subunit alpha in activated T-cells and allergic inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD124/IL-4R α mAb
Remarks/Alias	CD124; IL4RA; IL-4RA
Species	Human
GeneID (Human)	3566
GeneID	3566
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-232 of human CD124/IL-4R α (NP_000409.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P24394
Protein Weight	26kDa/88kDa/90kDa
Shipping	Ice bag